


Advanced laboratory class 2

FP2 - Nonlinear Optics - Second Harmonic Generation

Marco CANTERI
marco.canteri@student.uibk.ac.at

Innsbruck, November 17, 2017



Abstract

In this work we generated ultraviolet light at around 317 nm from a laser beam of 633 nm exploiting Second Harmonic Generation (SHG), a second order non linear effect of a potassium dihydrogen phosphate (KDP) crystal. We measured the power of the red laser as a function of the angle of a polarizer, then we studied the efficiency of the SHG with respect to the crystal angles and input power.

1 Theory

Inside a medium the relation between the polarization and the electric field E is in first approximation linear. If the intensity of the electric field is strong enough, the relation between the field and the polarization is no longer linear and we have to consider higher-order effects. In general we find

$$P_i = \varepsilon_0 \left(\sum_j \chi_{ij}^{(1)} E_j + \sum_{j,k} x_{ijk}^{(2)} E_j E_k + \sum_{j,k,l} \chi_{ijkl}^{(3)} E_j E_k E_l + \dots \right) \quad (1)$$

where the $\chi^{(i)}$ are tensors of rank $i + 1$ and represent the i -order susceptibility. This leads to an entire new class of phenomena. In fact it is possible, as we will show below, to excite a new electric field of frequency ω_2 with an electric field of frequency ω_1 .

From the Maxwell equations we can obtain [1] the following differential equation for the electric field

$$\left(\nabla^2 - \frac{1}{c^2} \frac{\partial^2}{\partial t^2} \right) E = \mu_0 \frac{\partial^2 P_{NL}}{\partial t^2} \quad (2)$$

where P_{NL} are the non linear terms of the polarization. We can see from this equation that the non linear polarization acts like a source for the electric field. Therefore, the electric field will oscillate at the same frequency of the polarization. Now, if we consider an exciting wave of frequency ω :

$$E = \frac{1}{2} (A(\omega) e^{i\omega t} + \text{c.c.}),$$

and we neglect any non linear order greater than two, we can see that P_{NL} , which goes with the square of the field, contains several terms of different frequency. We will focus on that which has a frequency of 2ω . Hence, we will look at generated light with frequency 2ω , this process is called Second Harmonic Generation.

In this experiment we used a KDP crystal to generate a frequency of 317 nm from a red laser of 633 nm. The main purpose was to study the efficiency of such process. The power of the generated light can be written [1] as

$$I_{2\omega} \propto \text{sinc}^2 \left(L \frac{\Delta k}{2} \right),$$

where $\Delta k = \frac{4\pi}{\lambda} (n_{2\omega} - n_\omega)$ is called the phase matching relation and L is the length of the crystal. In order to get the maximum power, it must holds $\Delta k = 0$. Therefore, the refractive index at frequency ω must be equal to the refractive index of frequency 2ω . In a normal crystal this never occurs due to the dispersion of light, but we can exploit birefringence of the crystal, i.e. the crystal has different refractive indexes for different light polarization. For an extraordinary wave with angle θ , the refractive index is given by the following formula [1]

$$\frac{1}{n^2(\theta, \omega)} = \frac{\cos^2 \theta}{n_o(\omega)} + \frac{\sin^2 \theta}{n_e(\omega)},$$

where n_o is the ordinary refractive index and n_e is the extraordinary refractive index of the crystal. Therefore, if the incoming and the generated waves have different polarization, they will also travel with different refractive index that can be used in order to compensate dispersion. There are two different type of phase matching, when the incoming waves have ordinary polarization it is called Type I phase matching, otherwise, if the waves have perpendicular polarization, it is called Type II phase matching. This condition can be visualized with the help of the so called refractive index ellipsoid. In the figure we can notice that there are some points where the extraordinary refractive index with frequency ω is equal to the ordinary refractive index with frequency 2ω , that is where there is phase matching.




Figure 1: Refractive index ellipsoid, red line refers to frequency ω , blue line to 2ω

2 Experiment setup




Figure 2: Experiment setup. A red laser is pumped into a KDP crystal to generate ultraviolet light at 317 nm (shown in blue in this figure) detected with a photomultiplier. L_1 and L_2 are lenses

The experiment setup is depicted in figure 2 and it consists of a Helium-Neon laser which output a light of 633 nm, followed by an half wave plate to rotate the polarization and a polarizer. Then the light is reflected with two mirrors and sent through a potassium dihydrogen phosphate crystal. Finally the light is detected with a photomultiplier powered with 700 V after it is collimated and filtered. The photomultiplier was connected to an oscilloscope where we visualized and saved the data. The crystal is mounted on a rotating platform, such that it was possible to change the angle of the crystal orientation in order to study the efficiency of SHG. We had two degrees of freedom corresponding to two different angles θ and ϕ as shown in the photo below.




Figure 3: Angles of the crystal on which we were able to act

3 Measurements and analysis

Before performing the measurements on SHG, we first measured the power of the laser for different angles of the polarizer. Measurements were taken with a powermeter, we took 10 seconds of measurements and we used the built-in function for the average value and the standard deviation. We started from angle 0° to 360° with a step of 10 degrees. As can be seen from the plot in figure 4, the data are in agreement with the theoretical law $P = P_{max} \cos^2(\theta)$, a fit has been done and it leads to $P_{max} = 7.24 \pm 0.07$ mW.




Figure 4: Power of the Helium Neon laser as a function of the polarizer angle. The error of the angle is due to the resolution of the polarizer, anyway the error bars are too small to be seen

Then, we set the polarizer at 180 degrees, so the light that went through the crystal was vertically polarized. Thus we measured the power of the second harmonic generation process with a photomultiplier. We chose a fixed angle $\theta = 11^\circ$, while we changed the angle ϕ . We took measurements from 0° to 360° with a step of 10 degrees. In the oscilloscope we acquired data with a scale of $500 \mu s$, we averaged these data and calculated the standard deviation in order to

get an error, while the error in the angle is only due to the resolution. The plot is shown in figure 5. We can notice two dips, one at around 130° and the other one at around 290° . The minimum in voltage corresponds to maximum in the power of the SHG. Therefore on these angles phase matching is achieved and we can observe the frequency doubling.




Figure 5: Photomultiplier signal of the second harmonic as a function of ϕ , experimental data are shown in point, while the line is only for eyes helping. The errors are too small and cannot be seen

After this measure, we kept ϕ fixed at $\phi = 290^\circ$ and we changed θ . We went from 37° to -2° due to limitation of the setup. The crystal and the lenses were too close to allow a large rotation of the crystal. Errors are calculated as before and we used the same scale on the oscilloscope.




Figure 6: Photomultiplier signal of the second harmonic as a function of θ , experimental data are shown in point, while the line is only for eyes helping.

The data are shown in figure 6, we can notice a dip around 15° , that is a maximum in the power of the second harmonic, i.e. the angle where phase matching is achieved. Finally we measured the power of the second harmonic as a function of the polarizer angle, i.e. the power of the He-Ne laser as shown in this analysis. We fixed $\phi = 290^\circ$ and $\theta = 11^\circ$, then we took the measurements with the same configurations of before and we evaluated the errors in the same way. The results are in figure 7. There is a peak around 180° as expected, in fact at this angle the incident light has the maximum power as can be seen in figure 4.




Figure 7: Photomultiplier signal of the second harmonic as a function of polarizer angle, experimental data are shown in point, while the line is only for eyes helping.

4 Summary and conclusion

In this experiment we used a Helium Neon laser to excite SHG inside a non linear crystal. First we checked the power of the source as a function of the polarizer angle and we found the expected cosine dependence. Then we studied the efficiency of the SGH process measuring its power as a function of the crystal angles. We obtained peaks in the power when the phase matching condition was achieved. This achievement is possible due to the birefringence of the crystal.

References

- [1] BAHAA E. A. SALEH, MALVIN CARL TEICH, *Fundamentals of photonics*, Wiley series in pure and applied optics, 1991, 1st edition
- [2] Fortgeschrittenenpraktikum 2, *Experiment FP2-07: Nonlinear Optics - Second Harmonic Generation*. SLAVA M. TZANOVA, KLEMENS SCHUPPERT.